



DASH® SARS-CoV-2 & Flu A/B Test External Control Log

Site Name: _____

Device Name: DASH® SARS-CoV-2 & Flu A/B Test _____

Test external controls once with each new shipment received and once for each untrained operator. Further controls may be tested local, state and/or federal regulations, accrediting groups, or your lab's standard Quality Control procedure.

Date	Operator's Initials	Test Lot #	Control Type Pos/Neg		Contol Lot/Exp	Results	Comments
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			
			Pos	Neg			

Reviewed by: _____

Date: _____

For support issues, please contact Customer Support at 1-888-992-DASH (3274)
or email customersupport@nuclein.com